

Aratohu Poto

CareSens™ NPOP

Pūnaha Aroturuki Huka Toto



1. Pūine Huka Toto CareSens N POP

Mata Whakaatu
Ka whakaatu huanga, karere

Pātene S
Ka whakakā/whakaweto rānei i te pūine, ka whakaū i ngā tipakonga tahua, ka huri mōhiohio hoki

Tauranga Ngaku Whakamātau
Kōkuhu te ngaku whakamātau ki konei

Pātene >>
Ka tipako mōhiohio, ka huri mōhiohio rānei

Tauranga Raraunga
Ka whakamahia hei whakawhiti raraunga i te pūine ki tētahi rorohiko whai taura

2. Ngaku Whakamātau CareSens N POP

Ngā pae whakapā
Āta pēhia te ngaku whakamātau, me ōna pae whakapā e anga ana ki runga, ki te tauranga ngaku whakamātau o te pūine

Pihanga whakaū
Titiro ki konei kia kite mēnā he rawaka te taura toto kua hoatu

Te tapa hei hoatu i te taura toto
Hoatu ki konei te taura toto hei whakamātau

3. Pūrere Wero

Pito Urutaua

Puringa Wero

Pūnaha Tuwha

Wero

Pihanga Whakaū Utanga

Pātene Tuku

Kāho Reti

4. Te whakarite i te Pūrere Wero

- Whakawiria ka tango i te pito urutaua.
- Āta kōkūhuna he wero hou ki te puringa wero. Pupuri paitia te wero. Āta whakawiria hei tango i te taupoki papare. Tiakina te taupoki hei poki anō i te wero i muri ake. Whakahokia te wero ki tōna taupoki.
- Tipakohia he hōhonu mai i te kotahi ki te whitu (1-7) ka hiahiaia i te pito urutaua o te pūrere wero. Kōwhiria he hōhonu mā te huri i te wāhi whakarunga o te pito urutaua kia hāngai te tau ki te pere. Ka taunakitia he tautuhinga tīmata o te toru (3).
- Whakakeua te pūrere wero mā te kume i te kāho reti kia ngoto te pūrere. Ka huri pūwhero te pihanga whakaū utanga.

5. Te whakarite i te Pūine me te Ngaku Whakamātau

- Horoia ōu ringa me te wāhi taurira. Opeopea ka āta whakamaroke.
- Kōkūhuna te ngaku whakamātau ko ngā pae whakapā e anga ana ki runga. Ka puta te tohu ▲ ki te mata.

6. Te Hoatu Toto

- Tikina he taura toto ki te pūrere wero.
- Putā rawa atu te tohu ▲ ki te mata, hoatu te taura toto ki te pito whāiti o te ngaku whakamātau kia pihō rā anō te pūine. Ki te mea kāore i kī te pihanga whakaū, ka puta he karere Er4.
- Ka puta te huanga whakamātau i muri i te tauanga iho o te pūine mai i te 5 ki te 1. Ka tiakina aunoatia te huanga ki te pūmahara o te pūine.

7. Te Maka i ngā Wero kua Whakamahia

- Pana whakamua te pūpei wero mā tō kōnui kia makaia ki tētahi ipu para koirora pūmate te wero kua oti te whakamahi.

8. Te Mata me te Karere

- Er1** Kua kuhuna he ngaku kua oti te whakamahi.
> Tukuruatia te whakamātau ki tētahi ngaku hōu.
- Er2** I hoatu he taura toto, taurira whakamātau rānei i mua i te putanga mai o te tohu ▲.
> Tukuruatia te whakamātau ki tētahi ngaku hōu ka tatari kia puta mai te tohu ▲ i mua i te hoatu i te taurira toto, te taurira whakamātau rānei.
- Er4** He piapiatanga nui rawa tō te taurira toto, he poto rānei te rōrahi.
> Tukuruatia te whakamātau i muri i te kuhunga o te ngaku hōu.
- Er5** Ka puta pea tēnei karere hapa ina whakamahia he ngaku whakamātau huka toto hē, kāhore ia te ngaku whakamātau huka toto CareSens N.
> Tukuruatia te whakamātau ki tētahi ngaku CareSens N.
- Er6** He raru kei te pūine.
> Kauga te pūine e whakamahia. Whakapā ki tō māngai tauhokohoko whaimana o i-SENS.
- HIC**
- LoL**

I te whakamātau i runga kē te pāmahana i te awhe mahinga.

I te whakamātau i raro kē te pāmahana i te awhe mahinga.

> Nuku ki tētahi wāhi kei roto te pāmahana i te awhe mahinga (5-50°C) ka tukuruia i te whakamātau ina tae te pūine me ngā ngaku whakamātau ki tētahi pāmahana i roto i te awhe mahinga.

Aratohu Poto

